



A. Model 4795



B. Model 595



C. Model 2010



D. Model 4416



E. Model 6045



F. Model 4086



G. Model 4074



H. Series 1427



Wheelchair Accessible

- A. Model 4795 Bariatric Hi-Lo Treatment Table with Power Backrest**  
Extra-wide 32" top accommodates patients up to 600 lbs.  
Power height motorized elevation from 20" wheelchair height to 29" high.  
Power backrest for positioning patient.
- B. Option 595 Digital 600 lb. Capacity Scale (Factory Installed)**  
Available option for Models 4795, 4074 and 4086. Pound/kilo display.
- C. Model 2010 Bariatric Footstool**  
Reinforced Steel cross-brace to accommodate up to 500 lbs.
- D. Model 4416 X-L Power-Back 600 lb. Capacity Exam Table**  
Super strong steel table for bariatric and general practice use.  
Power backrest for positioning patient.
- E. Model 6045 Bariatric 500 lb. Capacity Electric Tilt Table**  
Electric tilt mechanism with heavy-duty motor for up to 500 lb. patients.
- F. Model 4086 Bariatric 600 lb. Capacity Treatment Table**  
H-Brace with full width laminate shelf.
- G. Model 4074 Extra-Wide 800 lb. Capacity Bariatric Treatment Table**  
Extra-wide 34" upholstered top on cabinet style base.
- H. Series 1427 Bariatric Electric Mat Platform**  
Lift system can easily lift up to 750 lbs from 20" wheelchair height to 30" high.

Refer to our catalog or visit us on the web ([www.hausmann.com](http://www.hausmann.com)) for complete product features.