

LED ILLUMINATING EXAM LIGHT



Model 2179

*Ideal For General Examination**

Model 2179

LED Illuminating Exam Light

- Excellent Energy Savings Compared to Incandescent or Halogen Bulbs
- Green & Environmentally Friendly - No Lead and No Mercury
- Long Life – Up to 50,000 Hours
- Natural Lighting - Color Temperatures Up To 4,500 °K
- Super bright light, without heat
- No UV (Ultraviolet Light)
- Light Weight and Shock Resistant

*NOT INTENDED FOR PEDIATRIC, SURGERY OR DENTAL USE.

Technical Data:

Color Temperature:	4,500 °K
Color Rendition Index (CRI):	75 - 85
Illumination (Intensity):	≥ 12,000 lux (≥ 1,115 foot-candles)
Optimum Working Distance:	12" to 20"
Spot Size at Working Distance:	3" to 5" Ø
Type of Bulb:	Multiple Super Bright LED 3V/0.5W
Number of Bulbs:	8 LED Array
LED Est. Life Span:	50,000 hours
Electric Specifications:	Input: 120 VAC / 60 Hz
Power Consumption:	≤ 3.5 watts
Control:	ON/OFF Switch (No Dimmer)
Fixture O/A Height:	40" – 60"
Flexible Gooseneck Arm Length:	20"
LED bulbs are protected by a clear cover	
Fixture includes 5 castor base as Standard	