



# HAUSMANN BARIATRIC PRODUCTS

REHAB & THERAPY PRODUCT GROUPING

*BARIATRIC line designed specifically for obese patients.....*



750  
LBS

A. Model 1564



B. Model  
6045-709

500  
LBS



C. Model 6185

500  
LBS




1100  
LBS

D. Model 1429-57



400  
LBS

E. Model 4079



700  
LBS

F. Model 1440-57-113




500  
LBS

G. Model 6065



H. Model 5085

500  
LBS



600  
LBS

I. Model 1386

- A. Model 1564 Bariatric Ada Two-Sided Staircase:  
6"H steps on one side, 4"H steps on the other.
- B. Model 6045-709 Bariatric Electric Tilt Table:  
Electric tilt mechanism with heavy-duty motor. Extra-wide top.
- C. Model 6185 Bariatric 500 lb. Electric Stand-In Table:  
Power lifts patient from sitting to standing position.
- D. Series 1429 Heavy-Duty Electric Mat Platform:  
Motorized power height adjustment.
- E. Model 4079 Modality Treatment Table:  
Mid-section of table has full extension equipment compartment.
- F. Series 1440 Dual Lift Powermatic® Mat Platforms:  
Power height adjustment from 19" wheelchair height to 28"H.
- G. Model 6065 Electric Hi-Lo Bo-Bath Table:  
Combination treatment table and mat table.
- H. Model 5085 Bariatric Wood Transfer Board:  
5/8" thick with tapered ends.
- I. Model 1386 Bariatric Parallel Bars, Height & Width Adjustable:  
38" standard clearance between uprights.

CONTACT: →