




LONG TERM CARE PRODUCTS

MAT PLATFORMS



700 LBS A. Model 1440-57-113



600 LBS B. Model 1441-47-113



750 LBS C. Model 1427-57-121

- A. Series 1440 Dual Lift Powermatic® Mat Platforms** Power height adjustment from 19" wheelchair access to 28"H. (4 sizes)
 - B. Series 1441 Dual Lift Powermatic® Mat Platforms** Power height adjustment from 20" wheelchair access to 29"H. (3 sizes)
 - C. Series 1427 Bariatric Electric Mat Platforms** Extremely stable lift system easily lifts 750 lbs from 20" wheelchair height to 30"H. (4 sizes)
- Crank Hydraulic Mat Platform (Series 1425) also available (not shown). (4 Sizes)**

PARALLEL BARS



400 LBS D. Model 1393




400 LBS F. Model 1398

PULLEYS



G. Model 2629



600 LBS E. Model 1386

- D. Model 1393 Patented Ambulation & Mobility Platform** Height and width adjustable.
- E. Model 1386 Bariatric Parallel Bars, Height & Width Adjustable**
- F. Model 1398 Patented Mobility Platform with Electric Height Bars** Equipped with grommets holes for Ladder Steps and Balance Beam.
- Model 1396 Power Height Parallel Bars (not shown)** Same as 1398 less grommets holes.
- G. Model 2629 Duplex Pulley Weights** For shoulder, elbow and chest exercises.

TRAINING STAIRCASES



350 LBS H. Model 1566



350 LBS I. Model 1580



350 LBS J. Model 1570

- H. Model 1566 "Mini" Staircase** Space Saver can fit in 8 ft. high ceiling. (2) 6" high steps and (3) 4" high steps. 30" wide steps.
- I. Model 1580 ADA Convertible "Space Saver" Staircase** Compact footprint takes up less space. Can be installed in straight line or 90° corner arrangement.
- J. Model 1570 Patented 3-in-1 Training Staircase** Unique training staircase allows therapist to train patient on (3) progressive step levels.

CONTACT: →

TILT TABLES



350
LBS

K. Model 6058

STAND-IN TABLES

L. Model 6180

300
LBS



M. Model 6175



300
LBS

K. Model 6058 Wheelchair Accessible Hi-Lo Tilt Table Electric tilt mechanism operated with hands-free foot control. Comfortable 76" x 26" cushioned top with (3) safety straps.

L. Model 6180 Hi-Lo Stand-In Table with Electric Patient Lift can power lift 300 lb. patient from sitting to standing position.

M. Model 6175 Hi-Lo EconoLine Stand-In Table with electric patient lift motorized patient lift can elevate a 300 lb. person from sitting to standing position.

WORK/ACTIVITY TABLES



200
LBS

N. Model 4335



250
LBS

O. Model 4320



200
LBS

P. Model 4375



150
LBS

Q. Model S-740

N. Model 4335 Round Hydraulic Work Table Hydraulic foot activated mechanism raises table height.

O. Model 4320 Crank Hydraulic Work Table Hand crank operated system provides 14" height range.

P. Model 4375 Electric Work Table Direct-drive electric lift mechanism elevates top from 27" to 43" high.

Q. Model S-740 Econo Work Table 48"L x 32"W x 26" - 34"H. Model S-745 Econo Work Table (not shown) 60"L x 36"W x 26" - 34"H.

WEIGHT RACKS



300
LBS

R. Model 5560



300
LBS

S. Model 5569



400
LBS

T. Model 5562



U. Model 1671

R. Model 5560 Econo Weight Rack Storage for weights, full-view mirror. **Model 5560-100 Weight Rack plus accessories.**

S. Model 5569 Multi-Purpose Weight/Storage Rack 2-Door storage cabinet. **Model 5569-100 Weight/Storage Rack plus accessories.**

T. Model 5562 Mega Rack Weight storage, work surface and storage drawers. **Model 5562-100 Mega Rack plus accessories.**

U. Model 1671 Portable Mirror Distortion-free with ANSI "shatter-stop" safety backing.

MOBILE THERAPY CABINET



250
LBS

V. Model 6690 Cubex™ Therapy System on Wheels Bolt-board, work surface, pegboard pad, activity board, and two-door cabinet.

STOOLS

250
LBS

W. Model 2153



X. Model 2018



1000
LBS

Y. Model 2038



W. Model 2153 Economy Air-Lift Stool Single lever height control.

X. Model 2018 Bariatric XL Safe Step with anti-tip feet.

Y. Model 2038 Bariatric XL Safe Step with Handrail with anti-tip feet.

THERAPY STATION



300
LBS

Z. Model 5556 Mobile Therapy Station Grab bar and 4" locking casters. **Model 5556-100 Mobile Therapy Station plus accessories.**

