



PROFESSIONAL STOOLS & SEATING

Air-Lift Stools

250 LBS

- 16" diameter seat, 3 1/2" thick foam.
- 20" - 27" height range.
- Durable 5-leg, black composite base meets BIFMA and ANSI standards.
- (5) 2" twin wheel black nylon casters.
- Choice of (18) standard vinyl colors.



Model 2113

Model 2123

Model 2133

Model 2143



Model 2153



Model 2157-S



400 LBS Model 2151



Model 2154



Model 2156-707

Economy Air-Lift Stools

- 16" diameter seat, 3 1/2" thick foam.
- 17"-22" or 20"-27" height range.
- (5) 2" twin wheel black nylon casters.

250 LBS

Heavy-Duty Air-Lift Stool

- 400 lbs. weight capacity.
- Pneumatic lift with "D" control handle.

Dual Foot & Hand Operated Stool

- 360° ring hand control & circular metal foot control.

250 LBS

Economy Air-Lift Stool

- Air-spring pneumatic height adjustment.
- 360° ring hand control.

250 LBS



5-Caster Revolving Stool Model 2150

250 LBS



Revolving Stool Model 2116

250 LBS



P.T. "Scooter" Stool Model 2300

250 LBS



A. 2165

750 LBS



B. 2158

250 LBS

A. Oak Bariatric Arm Chair

- 28"W x 19"D x 19"H seat

B. Arm Chair

- Steel frame with upholstered seat and backrest

C. Side Chair

- Steel frame with upholstered seat and backrest



C. 2159

250 LBS



Foot Stool with Safety Handrail Model 2030

300 LBS



Foot Stool Model 2000

Bariatric XL Safe Step with Handrail Model 2038

1000 LBS

Bariatric XL Safe Step Model 2018



1000 LBS

CONTACT: