



QUALITY-LINE TREATMENT TABLES



A. Model 4002

H-Brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.



B. Model 4024



C. Model 4041
Shown with option: 05 - Headrest

A. Series 4002

- High pressure laminate legs, apron, and stretchers.



MODEL	LENGTH	WIDTH	HEIGHT
4002*	72"	24"	31"
4002-027	72"	27"	31"
4002-030*	72"	30"	31"
4002-030-078	78"	30"	31"

B. Series 4024

- H-Brace with laminate storage shelf.
- High pressure laminate legs, apron, and stretchers.



MODEL	LENGTH	WIDTH	HEIGHT
4024*	72"	24"	31"
4024-027	72"	27"	31"
4024-030*	72"	30"	31"
4024-030-078	78"	30"	31"

C. Series 4041

- H-Brace with laminate shelf and open-shelf storage cabinet.
- High pressure laminate legs, apron, and stretchers.



MODEL	LENGTH	WIDTH	HEIGHT
4041	72"	24"	31"
4041-027	72"	27"	31"
4041-030	72"	30"	31"
4041-030-078	78"	30"	31"

D. Series 4043

- Fully-enclosed, all-laminate storage cabinet with sliding laminate doors.
- Interior has divider and (1) adjustable shelf.
- Black laminated plywood base.



MODEL	LENGTH	WIDTH	HEIGHT
4043	72"	24"	31"
4043-027	72"	27"	31"
4043-030*	72"	30"	31"
4043-030-078	78"	30"	31"



D. Model 4043

COMMON FEATURES

QUALITY LINE TREATMENT TABLES — SERIES 4002, 4024, 4041, 4043:



- 2" High-Density urethane foam top.
- Resist-All™ vinyl upholstery finish.
- Laminate: Natural Oak, Wild Cherry, or Folkstone Gray
- Upholstery – Choice of (18) standard vinyl colors

OPTIONS

FOR 4000 SERIES TABLES

- *00 – Special Table Height (N/A Model 4043)
- *04 – Table Safety Strap
- *05 – Headrest, adjustable
- *06 – Backrest, adjustable
- *07 – Nose Cutout
- *08 – Drawer (N/A with option 69)
- *15 – Safety Rail (only Model 4043)
- *19 – Slide Shelf, laminate
- 23 – Paper Dispenser
- 25 – Paper Cutter
- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small
- *69 – Gas Spring Backrest
- *99L – Left Hand Table Model

*FACTORY INSTALLED

ENHANCED LEG HARDWARE SYSTEM



Attach or detach legs quickly and easily with heavy-duty 5/16" bolts into embedded metal threaded inserts for durable connection. Easily tighten and re-tighten as needed.

CONTACT: