



a division of
Dynatronics

BARIATRIC MAT PLATFORM



MODEL 1427-57-121

Series 1427 Bariatric Mat Platforms



- Electrically activated hydraulic cylinders are built into each (4) corner steel leg.
 - Motorized lift system can securely lift 750 lbs. from 20" wheelchair height to 30" high.
 - "Hands-Free" foot control operation.
 - Seamless heavy duty tear-resistant nylon-reinforced vinyl.
 - 2" High-Density urethane foam top.
 - Steel frame and 4 corner legs have durable black powder coat finish.
 - Legs have levelers for uneven floors.
 - 4-Leg "Open" design is patient sling-lift accessible on all four sides.
- * You must specify upholstery color with order.

UPHOLSTERY COLORS

113 Blue

119 Gray

SPECIAL ORDER COLORS

121 Green

122 Black

MODEL	WIDTH	LENGTH	HEIGHT
1427-46	4'	6'	20" - 30"
1427-47	4'	7'	20" - 30"
1427-57	5'	7'	20" - 30"
1427-68	6'	8'	20" - 30"