




Series 1440

Dual Lift Powermatic® Mat Platforms

- Power height adjustment from 19" wheelchair height to 28"H.
- Two direct drive, synchronized lift mechanisms for unsurpassed power.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" High-Density urethane foam padding for comfort and rounded corners for safety.
- Black laminate base with satin aluminum trim and levelers.
- "Hands-Free" foot control.
- Hospital grade plug. 110V, 60Hz, 2.8 Amps. U.L. Listed.
- You must specify upholstery color with order.



MODEL 1440-57-113

MODEL	WIDTH	LENGTH	HEIGHT
1440-46	4'	6'	19" - 28"
1440-47	4'	7'	19" - 28"
1440-57	5'	7'	19" - 28"
1440-68	6'	8'	19" - 28"



MODEL 1440-47-119



UPHOLSTERY COLORS



113 Blue



119 Gray

SPECIAL ORDER COLORS (NOT QUICKSHIP)



121 Green



122 Black

