



# HEAVY DUTY COLOR-CODED NESTED FOOTSTOOLS



MODEL 2418

## Model 2418 Heavy Duty Color-Coded Nested Footstools



- 3½", 6", and 8½" color coded heights.
- 19" x 12" platforms with Anti-Slip tread.
- Heavy duty 1,000 lb. capacity.
- Light weight construction – Each stool is ONLY 5 lbs.
- All-aluminum welded construction.
- Wide Anti-Tip bases with rubber feet.
- "Outrigger" legs are very stable – Great for geriatric & pediatric use!

MODEL	LENGTH	WIDTH	HEIGHT	PLATFORM L x W
2418 (NESTED SET)	20"	27½"	8½"	
Yellow	17"	25½"	3½"	19"x 12"
Red	17"	26½"	6"	19"x 12"
Blue	17"	27½"	8½"	19"x 12"