



a division of  
**Dynatronics**

# EASY TO LIFT OVER BED TABLE



**MODEL 3415**

## Model 3415 Easy To Lift Over Bed Table



- All metal design with 1" laminate top for extra strength and durability.
- Ultra-smooth gas lift height adjustment locks table at desired height.
- Height adjustable 30" - 43".
- Durable and user-friendly functionality for all healthcare facilities.
- 32" x 15" x 1" Oak laminate top.
- Fits standard hospital beds and nursing home beds.
- "C" base design to accommodate both wheelchair and chair access.
- Elegant design and finished metal for easy cleaning and maintenance.
- Four 2" locking casters for easy roll and security.
- Base height requires only 3" clearance.

32" x 15" x 1" Oak laminate top



MODEL	LENGTH	WIDTH	HEIGHT
3415	32"	15"	30"- 43"
	81cm	38cm	76-109cm