



# TORQUE THERAPY TABLE

MODEL 4286

*Your extra pair of hands.*




## STABILIZE




Iliotibial Band (ITB) Stretch

## ISOLATE



Paraspinal Roll

## MANIPULATE



Lumbar Traction

## FIXATE



Cervical Traction

### Model 4286 Torque Therapy Table



- Table is designed to aid in stretching, muscle isolation, joint distraction, manual spinal manipulation, and orthopedic special testing.
- Twelve strategically placed holes allow for placement of two different size movable bolsters.
- Includes two 6" round x 12" tall bolsters and two 3" round x 12" tall bolsters designed to be placed vertically on the table.
- Six legs and plywood top for added strength and 600 lb weight capacity.
- Natural Oak high pressure laminate legs, aprons, and H-brace.
- 1.5" reversible blue vinyl coated high-density mat measures 30"W x 78"L.
- Full length shelf for storage of bolsters when not in use.
- Legs have protective black vinyl scuff caps.
- Patent pending design.



SCAN FOR DEMO VIDEO

