



HAUSMANN POWER MAT PLATFORMS

BARIATRIC HEAVY-DUTY MAT PLATFORM

lifts 1100 lbs. from 20" wheelchair height to 30" high



1429-68-113

electrically activated hydraulic cylinders are built into each corner leg

4 leg "open" design is patient sling-lift accessible on all four sides



Available in two sizes:

MODEL	WIDTH	LENGTH	HEIGHT
1429-57	5'	7'	20"-30"
1429-68	6'	8'	20"-30"




BARIATRIC MAT PLATFORM



1427-57-121

DUAL LIFT POWERMATIC® MAT PLATFORM with Black Laminate Top



1449-57

shown with option 3002-57 mat


ECONO ELECTRIC HI-LO MAT PLATFORM



1441-47-113

DUAL LIFT POWERMATIC® MAT PLATFORM

two direct-drive, synchronized lift mechanisms for unsurpassed power



adjusts from 19" wheelchair height to 28" h.




1440-57-113




hands-free foot control


MAT NYLON-REINFORCED VINYL COLORS:



113 BLUE



119 GRAY



121 GREEN



122 BLACK

TURN OVER
FOR MORE MODELS



CONTACT:



HAUSMANN MAT PLATFORMS

CRANK HYDRAULIC MAT PLATFORM

Shown with optional Bolsters #41 and #43

2" high-density urethane foam
top upholstered with seamless
tear-resistant reinforced vinyl

hydraulic
cylinders are
built into each
corner leg

black powder
coated steel
frame

1425-68

easy-to-use
crank hydraulic
system can lift
750 lbs. from 20"
wheelchair height
to 30" high

4-Leg "open" design is
patient sling-lift accessible
on all four sides

leg levelers

TURN OVER
FOR MORE MODELS

Available in four sizes:

MODEL	WIDTH	LENGTH	HEIGHT
1425-46	4'	6'	20"-30"
1425-47	4'	7'	20"-30"
1425-57	5'	7'	20"-30"
1425-68	6'	8'	20"-30"

750
LBS

UPHOLSTERED MAT PLATFORM



1412-47

BARIATRIC SPACE-SAVER THERAPY PLATFORM



1418

ECONO MAT PLATFORM



750
LBS

S-205

SPACE-SAVER MAT PLATFORM

upholstered with seamless heavy duty,
tear-resistant nylon-reinforced vinyl

fold away
design saves
floor space



three folding high
pressure laminate legs

1420

18" height is convenient for
therapist and wheelchair
accessible for patient

sturdy,
floor-supported
mounting plate

BARIATRIC WALL FOLDING THERAPY PLATFORM



1000
LBS

1422

Visit our website or refer to current catalog pages 40-45 for full product details and available top sizes: www.Hausmann.com