



# SILVER STAR TABLES & EQUIPMENT


## A. Silver Star Tables

- Oak laminate on all surfaces.
- Solid end panels with full-width storage shelf.
- Upholstered flat top has 2" urethane foam and choice of 706 Blue or 709 Gray vinyl.
- Ships ready-to-assemble.
- 325 lbs. weight capacity.



A. Model S-110

MODEL	LENGTH	WIDTH	HEIGHT	UPH.COLOR
S-109	72"	27"	31"	706 Blue
S-110	72"	27"	31"	709 Gray



### Options FOR S-109, S-110

- 27 – Paper Dispenser & Cutter Combo
- 32 – Pillow, full size
- 35 – Pillow, small




B. Model S-205

## B. Silver Star Mat Platforms



- Sturdy panel legs with wear-resistant Folkstone Gray laminate
- 2" foam standard, 18" height.
- Upholstered in seamless nylon reinforced 113 Blue or 119 Gray.

MODEL	COLOR	WIDTH	LENGTH	HEIGHT
S-201	Blue	4'	7'	18"
S-202	Gray	4'	7'	18"
S-205	Blue	5'	7'	18"
S-206	Gray	5'	7'	18"



C. Model S-601

## C. Silver Star Two-Sided Staircase



- Two 6" steps and three 4" steps.
- Steps and platform have "anti-slip" safety treads.
- Platform is 30" W x 24" D x 12" H.
- Compact size can fit in 8' ceilings.
- Folkstone Gray laminate finish.
- Ergonomic hardwood handrails.

MODEL	LENGTH	WIDTH	HEIGHT
S-601	52"	30"	12"

CONTACT: →



D. Model S-410



E. Model S-420-100

### D. Silver Star Econo Mirror

- Full-view glass mirror with ANSI safety backing.
- Gray laminate with (4) 2 1/2" casters.

MODEL	WIDTH	DEPTH	HEIGHT
S-410	24"	18"	70"

### E. S-420 Silver Star Mirror Rack

- Same as S-410 but with Pegboard back and (11) hooks.

MODEL	WIDTH	DEPTH	HEIGHT
S-420	24"	18"	70"

S-420-100 Consists of: (1) S-420 Mirror Rack,  
(1) 5581 Set of (16) Huger Weights

### F. Silver Star Wall Mounted Overhead Pulley


- For light duty upper extremity exercise.
- 10" x 8" gray laminate wall-mounting plate.
- Equipped with (2) swivel rollers, nylon rope and two handles.



F. Model S-950

MODEL	DEPTH	WIDTH	HEIGHT
S-950	10"	16"	8"

### G. Silver Star Platform Mounted Parallel Bars

- 1 1/2" diameter one-piece stainless steel handrails for easy use and durability. 
- Heavy gauge steel base plates, uprights and fittings. Black powder coated chip and wear resistant finish.
- Each upright telescopes up in 1/2" increments and locks into (10) height positions with fail-safe ball-tip locking pin.
- Folkstone Gray laminate platform with tapered ends and anti-slip rubber treads.
- Height and width adjustable bars.

 350 LBS



G. Model S-320

MODEL	LENGTH	WIDTH*	HEIGHT
S-320	10'	15"- 28"	28"- 41"
S-321	7'	15"- 28"	28"- 41"

\*platform is 36" wide



H. Model S-740

### H. Silver Star Work Table

- Sturdy rectangular work table has no apron to accommodate wheelchair patient.
- Durable Gray high pressure laminate 1 1/8" laminate top with T-molding edge.
- (4) adjustable height metal legs.
- Hand control knob on each leg allows 1" adjustments from 26" to 34".

MODEL	LENGTH	WIDTH	HEIGHT
S-740	48"	32"	26"- 34"
S-745 (not shown)	60"	36"	26"- 34"