



# NEW STANDARD SKUS:

A GUIDE TO OUR MOST POPULAR PRODUCTS

Hausmann's new Standard SKUs\* were carefully chosen. These SKUs represent our most popular products, completely configured, including colors, sizes and options. These SKUs cannot be altered or substituted. Our goal is to "Make it Easy" for our customers to order the high-quality products they have been choosing for years.

This catalog only contains Standard SKUs. Customers can find Custom SKUs with available colors and options online at: [www.Hausmann.com/products](http://www.Hausmann.com/products)

## TABLE OF CONTENTS

3	Mat Platform
5	Treatment and Specialty Tables
13	Treatment and Power Tables
16	Pediatric Equipment
17	Recovery Couches
18	PROTEAM™ Multi-Seat Taping Stations
19	PROTEAM™ Treatment Tables
27	Wood Treatment Tables
28	HTT Tables

\*In this catalog, some Standard SKUs are not represented by our product images.

# MAT PLATFORM

## BARIATRIC UPHOLSTERED MAT PLATFORM



- All models have six (6) legs for increased strength and stability.
- 2" thick high-density urethane foam.
- Upholstered in seamless heavy duty, tear-resistant nylon-reinforced vinyl.

SKU	W	L	LAMINATE	VINYL
1412-46-L01V38	4'	6'	Natural Oak	Black
1412-47-L01V38	4'	7'	Natural Oak	Black
1412-57-L01V38	5'	7'	Natural Oak	Black
1412-68-L01V38	6'	8'	Natural Oak	Black

## SPACE SAVER WALL FOLDING MAT PLATFORM



- Fold-away design saves floor space.
- Sturdy, floor supported mounting plate.
- Three (3) folding high pressure laminate legs.



SKU	W	L	LAMINATE	VINYL
1420-47-L01V35	4'	7'	Natural Oak	Blue
1421-57-L01V38	5'	7'	Natural Oak	Black

# MAT PLATFORM

## BARIATRIC SPACE SAVER WALL FOLDING MAT PLATFORM



- Fold-away design saves floor space.
- 2" High-density urethane foam top.
- Sturdy, floor-supported wall mounting plate for added strength.

SKU	W	L	LAMINATE	VINYL
1422-3078-L01V20	30"	78"	Natural Oak	Dove Gray
1422-3078-L01V23	30"	78"	Natural Oak	Black

# TREATMENT AND SPECIALTY TABLES

## H-BRACE TREATMENT TABLE



- H-Brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.
- High-pressure laminate legs, apron and stretchers.
- 2" High-density urethane foam top.

SKU	W	L	LAMINATE	VINYL	OPTIONS
4002-2472-L01V26	24"	72"	Natural Oak	Royal Blue	
4002-2472-L02V17-A27	24"	72"	Folkstone Gray	Gray	Paper Dispenser & Cutter Combo
4002-3072-L02V17-A27	30"	72"	Folkstone Gray	Gray	Paper Dispenser & Cutter Combo
4002-3072-L02V19	30"	72"	Folkstone Gray	Navy Blue	
4002-3078-L01V16-A69	30"	78"	Natural Oak	Black	Gas Spring Backrest
4002-3078-L01V19	30"	78"	Natural Oak	Navy Blue	

# TREATMENT AND SPECIALTY TABLES

## H-BRACE TREATMENT TABLE WITH STORAGE SHELF



- H-Brace with laminate storage shelf.
- High-pressure laminate legs, apron and stretchers.
- 2" High-density foam top.

SKU	W	L	LAMINATE	VINYL	OPTIONS
4024-2472-L01V23-A27	24"	72"	Natural Oak	Black	Paper Dispenser & Cutter Combo
4024-2472-L02V15	24"	72"	Folkstone Gray	Slate Blue	
4024-2772-L01V12-A06A27	27"	72"	Natural Oak	Brown	Backrest, Paper Dispenser & Cutter Combo
4024-3072-L02V23	30"	72"	Folkstone Gray	Black	
4024-3078-L01V23-A06	30"	78"	Natural Oak	Black	Backrest

# TREATMENT AND SPECIALTY TABLES

## H-BRACE TREATMENT TABLE WITH STORAGE SHELF AND OPEN CABINET STORAGE



- H-Brace with laminate shelf and open-shelf cabinet.
- High-pressure laminate legs, apron and stretchers.
- 2" High-density foam top.

SKU	W	L	LAMINATE	VINYL	OPTIONS
4041-2472-L02V23	24"	72"	Folkstone Gray	Black	
4041-2472-L03V23	24"	72"	Wild Cherry	Black	
4041-2772-L01V17-A06A19	27"	72"	Natural Oak	Gray	Backrest, Sliding Shelf
4041-3072-L01V26-A27	30"	72"	Natural Oak	Royal Blue	Paper Dispenser & Cutter Combo
4041-3078-L01V15	30"	78"	Natural Oak	Slate Blue	

# TREATMENT AND SPECIALTY TABLES

## 24"x72" STORAGE TREATMENT TABLE WITH SLIDING CABINET DOORS



- Fully-enclosed, all-laminate storage cabinet with sliding laminate doors.
- Interior has divider and one (1) adjustable shelf.
- Black laminate plywood base.

SKU	W	L	LAMINATE	VINYL	OPTIONS
4043-2472-L03V15-A27	24"	72"	Wild Cherry	Slate Blue	Paper Dispenser & Cutter Combo
4043-2472-L04V17-A27	24"	72"	Maple	Gray	Paper Dispenser & Cutter Combo

## 34"x74" TWO PERSON CONVERTIBLE TAPING AND TREATMENT TABLE

**PROTEAM**  
BY HAUSMANN

- Solid one piece construction (not modular).
- Two (2) padded taping sections 26" wide x 34" deep x 36" high with one (1) removable middle top pad.
- Just add middle pad to convert taping table into a treatment table.



SKU	W	L	LAMINATE	VINYL
4061-3474-L01V19	34"	74"	Natural Oak	Navy Blue
4061-3474-L01V23	34"	74"	Natural Oak	Black



# TREATMENT AND SPECIALTY TABLES

## 27"x78" SINGLE PERSON CONVERTIBLE TAPING AND TREATMENT TABLE WITH GAS SPRING BACKREST



**PROTEAM**  
BY HAUSMANN

- Multipurpose table functions both as a taping bench with backrest and as a treatment table.
- To use as a taping bench, raise the gas spring padded backrest up to 75° with recessed control handles on each side of the table.
- Just lower the backrest and the table quickly converts into a full-length treatment table.

SKU	W	L	LAMINATE	VINYL
4062-2778-L01V26	27"	78"	Natural Oak	Royal Blue

## 30"x78" ORTHOPEDIC CASTING TABLE WITH STORAGE CABINETS



- Three (3) storage sections each with laminate interiors and one (1) adjustable shelf.
- 2" High-density foam top.

SKU	W	L	LAMINATE	VINYL
4073-3078-L02V20	30"	78"	Folkstone Gray	Dove Gray

# TREATMENT AND SPECIALTY TABLES

## 24"x72" SPACE SAVER WALL FOLDING CHANGING TABLE



- Fold-away design saves floor space.
- 2" High-density foam upholstered top.
- Sturdy, floor-supported wall mounting plate for added strength.

SKU	W	L	LAMINATE	VINYL
4075-2472-L01V15	24"	72"	Natural Oak	Slate Blue
4075-2472-L01V17	24"	72"	Natural Oak	Gray
4075-2472-L01V19	24"	72"	Natural Oak	Navy Blue

## 30"x78" TREATMENT TABLE WITH STORAGE CABINETS AND MODALITY DRAWER

- Designed to save time and effort with modality equipment at your fingertips.
- Spacious upholstered top with 2" thick high-density foam.
- Two (2) fully enclosed door cabinets with adjustable shelf.



SKU	W	L	LAMINATE	VINYL
4079-3078-L01V23	30"	78"	Natural Oak	Black
4079-3078-L01V26	30"	78"	Natural Oak	Royal Blue

# TREATMENT AND SPECIALTY TABLES

## 30"x78" H-BRACE BARIATRIC TREATMENT TABLE WITH STORAGE SHELF



- Six (6) high-pressure laminate legs and a sturdy plywood sub-top.
- H-Brace with full-width laminate shelf. 4" clearance under shelf.
- 2" High-density urethane foam top.

SKU	W	L	LAMINATE	VINYL
4086-3078-L01V15	30"	78"	Natural Oak	Slate Blue
4086-3078-L01V17	30"	78"	Natural Oak	Gray

## 24"x52-72" H-BRACE SPACE-SAVER TREATMENT TABLE

- Hausmann Space-Saver tables meet your needs for a full-size table in a compact space. When closed, tables are only 52" long. When open, they extend to a conventional 72".
- Full-width 20" long foot end folds down to reduce table length to 52".
- Positive action locking device allows the foot end to be easily elevated into three positions: 55°, 75° and 90°.



SKU	W	L	LAMINATE	VINYL
4102-2457-L02V17	24"	72"	Folkstone Gray	Gray
4102-2472-L03V26	24"	72"	Wild Cherry	Royal Blue

# TREATMENT AND SPECIALTY TABLES

## 24"x52-72" H-BRACE SPACE-SAVER TREATMENT TABLE WITH STORAGE SHELF



- Hausmann Space-Saver Tables meet your needs for a full-size table in a compact space. When closed, tables are only 52" long. When open, they extend to a conventional 72".
- Full-width 20" long foot end folds down to reduce table length to 52".
- Positive action locking device allows the foot end to be easily elevated into three positions: 55°, 75° and 90°.

SKU	W	L	LAMINATE	VINYL	OPTIONS
4124-2472-L01V17-A06	24"	72"	Natural Oak	Gray	Backrest
4124-2472-L03V20-A06A27	24"	72"	Wild Cherry	Dove Gray	Backrest, Paper Dispenser & Cutter Combo

## 24"x52-72" SPACE-SAVER TREATMENT TABLE WITH STORAGE CABINET

- Space-Saver Tables meet your needs for a full-size table in a compact space. When closed, tables are only 52" long. When open, they extend to a conventional 72".
- Full-width 20" long foot end folds down to reduce table length to 52".
- Positive action locking device allows the foot end to be easily elevated into three positions: 55°, 75° and 90°.



SKU	W	L	LAMINATE	VINYL	OPTIONS
4143-2472-L01V12-A05A27	24"	72"	Natural Oak	Brown	Paper Dispenser & Cutter Combo, Headrest



# TREATMENT AND POWER TABLES

## MAMMOTH 2: 2-SECTION HI-LO EXAMINATION TABLE



- Tables take up less operational space with no end-to-end movement
- Heavy duty independently lockable casters
- 2.5" tabletop featuring rounded corners and firm foam padding
- 550-pound weight capacity

SKU	W	L	METAL	VINYL
7702-L00V20	30"	76"	Cool Gray	Dove Gray

## MAMMOTH 3: 3-SECTION HI-LO TREATMENT TABLE

- Tables take up less operational space with no end-to-end movement
- Heavy duty independently lockable casters
- 2.5" tabletop featuring rounded corners and firm foam padding
- 550-pound weight capacity



SKU	W	L	METAL	VINYL
7703-L00V23	30"	76"	Cool Gray	Black
7703-L00V20	30"	76"	Cool Gray	Dove Gray
7703-L00V19	30"	76"	Cool Gray	Navy Blue

# TREATMENT AND POWER TABLES

## 30"x72" HANDS FREE ELECTRIC HI-LO TREATMENT TABLE



- Adjustable height motorized hydraulic cylinders are built into each four (4) corner steel leg.
- "Hands Free" foot control activation.
- Apron and shelf are almond laminate.

SKU	W	L	METAL	VINYL	OPTIONS
4707-3072-L00V19-A27A06	30"	72"	White	Navy Blue	Backrest, Paper Dispenser & Cutter Combo

## 30"x70" ALL-PURPOSE STORAGE TREATMENT TABLE WITH 3 DRAWERS AND CABINETS

- Three (3) all laminate drawers on nylon rollers.
- Three (3) fully-enclosed door cabinets.
- 2" High-density urethane foam top.



SKU	W	L	LAMINATE	VINYL	OPTIONS
4835-3078-L01V23-A06	30"	78"	Natural Oak	Black	Backrest
4835-3078-L02V15	30"	78"	Folkstone	Slate Blue	

# TREATMENT AND POWER TABLES

## 24"x41-72" EXTENDED 2-IN-1 TREATMENT AND EXAM TABLE WITH STORAGE CABINETS



- Hausmann 2-in-1 Tables save space and costs by providing a readily adjustable examination and treatment table in a single, attractive and functional unit.
- Extended, the table is 72" long with a 14" adjustable headrest, a 17" adjustable footrest and a 41" center section.
- Adjustable length stirrups are included as standard equipment.

SKU	W	L	LAMINATE	VINYL	OPTIONS
4643-2472-L02V13-A27	24"	72"	Folkstone Gray	Rodeo Tan	Paper Dispenser & Cutter Combo

## PEDIATRIC EQUIPMENT

### 24"x43" ECONO-LINE™ PEDIATRIC DIGITAL SCALE WITH DRAWERS AND STORAGE SHELF



- Digital scale with lb-oz or kg readout.
- Reads in 1 oz. increments.

SKU	W	L	LAMINATE	VINYL
4941-2443-L02V41	24"	43"	Folkstone Gray	Flint Gray

### 24"x43" ECONO-LINE™ PEDIATRIC DIGITAL SCALE WITH STORAGE CABINETS AND DRAWERS



- Digital scale with lb-oz or kg readout.
- Reads in 1 oz increments.

SKU	W	L	LAMINATE	VINYL
4945-2443-L02V26	24"	43"	Folkstone Gray	Royal Blue
4945-2443-L03V40	24"	43"	Wild Cherry	Indigo Blue



# RECOVERY COUCHES

## 24"x72" ECONO-LINE™ RECOVERY COUCH WITH SLIDING CABINET DRAWERS



- Wood has no added urea-formaldehyde for safer indoor air quality.
- CARB II compliant.
- Nature-friendly urethane.

SKU	W	L	LAMINATE	VINYL	OPTIONS
7031-2472-L02V40-A05A27	24"	72"	Folkstone Gray	Indigo Blue	Headrest, Paper Dispenser & Cutter Combo

# PROTEAM™ MULTI-SEAT TAPING STATIONS

## MODULAR TAPING STATIONS



**PROTEAM**  
BY HAUSMANN

- Upholstered Taping Modules are 36" or 42"H.
- Storage Modules with laminate tops are 4" lower.
- Durable dowel and glue laminated wood construction for long term structural integrity.
- Self-closing laminate drawer with nylon rollers.
- Doors have concealed, adjustable, self-closing 110° Euro-hinges.

SKU	W	L	LAMINATE	VINYL	OPTIONS
A521-4236-L01V19-A26	36"	72"	Natural Oak	Navy Blue	(3) Fixed 75° Backrest
A521-4236-L01V23-A26A24A18A23E	36"	72"	Natural Oak	Black	(2) Fixed 75° Backrest, (1) Waste Hamper, (1) 3 Drawers, (2) Full Extension Pull-Out, (2) Emboss Logo
A531-4236-L01V19-A26	36"	120"	Natural Oak	Navy Blue	(3) Fixed 75° Backrest
A531-4236-L01V23-A26A24A18A23E	36"	120"	Natural Oak	Black	(3) Fixed 75° Backrest, (2) Waste Hamper, (2) 3 Drawers, (3) Full Extension Pull-Out, (3) Emboss Logo
A541-4236-L01V19-A26	36"	168"	Natural Oak	Navy Blue	(3) Fixed 75° Backrest
A541-4236-L01V23-A26A24A18A23E	36"	168"	Natural Oak	Black	(4) Fixed 75° Backrest, (3) Waste Hamper, (3) 3 Drawers, (4) Full Extension Pull-Out, (4) Emboss Logo

# PROTEAM™ TREATMENT TABLES

## 30"x78" GAS SPRING LEG REST RECOVERY TABLE WITH SLIDE PULLOUT DRAWERS



**PROTEAM**  
BY HAUSMANN

- Foot end of table has gas spring elevating leg rest (up to 75°) and a full extension modality drawer.
- Middle section of table has a side pull-out drawer with four (4) drying rods.
- Head section of table has two (2) hinged doors with interior adjustable shelf offering concealed storage.

SKU	W	L	LAMINATE	VINYL
A9021-3078-L03V23	30"	78"	Wild Cherry	Black

## 30"x78" H-BRACE TREATMENT TABLE WITH STORAGE SHELF

**PROTEAM**  
BY HAUSMANN

- 30" x 78" upholstered top with comfortable 2" thick high-density foam.
- High-pressure laminate legs, apron and H-Brace.
- Enhanced leg hardware system.



SKU	W	L	LAMINATE	VINYL
A9024-3078-L01V23	30"	78"	Natural Oak	Black

# PROTEAM™ TREATMENT TABLES

## 30"x78" TREATMENT TABLE WITH STORAGE SHELF AND OPEN CABINET STORAGE



**PROTEAM**  
BY HAUSMANN

- 30" x 78" upholstered top with comfortable 2" thick high-density foam.
- Laminate shelf and open storage cabinet.
- High-pressure laminate legs, apron and H-Brace.

SKU	W	L	LAMINATE	VINYL
A9041-3078-L01V19	30"	78"	Natural Oak	Navy Blue
A9041-3078-L01V23	30"	78"	Natural Oak	Black

## 30"x78" SPLIT LEG TREATMENT TABLE WITH OPEN CABINET STORAGE

**PROTEAM**  
BY HAUSMANN

- Open shelf storage, accessible from both sides of the table.
- Extremely versatile and functional split leg design.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.




SKU	W	L	LAMINATE	VINYL	OPTIONS
A9065-3078-L01V19-A69	30"	78"	Natural Oak	Navy	Gas Spring Backrest
A9065-3078-L02V23-A69	30"	78"	Natural Oak	Black	Gas Spring Backrest



# PROTEAM™ TREATMENT TABLES

## 30"x78" SPLIT LEG TREATMENT TABLE WITH STORAGE CABINETS



**PROTEAM**  
BY HAUSMANN

- Six (6) hinged doors with interior shelf offer concealed supply storage.
- Solid finished back.
- Extremely versatile and functional split leg design.

SKU	W	L	LAMINATE	VINYL	OPTIONS
A9068-3078-L01V19-A69	30"	78"	Natural Oak	Navy Blue	Gas Spring Backrest

## 30"x78" KNEE AND HIP FLEXION SPLIT LEG TREATMENT TABLE WITH STORAGE CABINETS

**PROTEAM**  
BY HAUSMANN

- Extremely versatile and functional split leg design.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Each padded leg rest can be adjusted with hip and knee joints locked into 45° angle with easy to operate hinge mechanism.



SKU	W	L	LAMINATE	VINYL
A9071-3078-L01V19	30"	78"	Natural Oak	Navy Blue

# PROTEAM™ TREATMENT TABLES

## 30"x78" TREATMENT TABLE WITH STORAGE CABINETS



**PROTEAM**  
BY HAUSMANN

- Three (3) two-door cabinets with laminate interiors and adjustable shelf.
- Solid finished back.
- Comfortable 2" thick high-density foam upholstered top.

SKU	W	L	LAMINATE	VINYL
A9073-3078-L01V23	30"	78"	Natural Oak	Black

# PROTEAM™ TREATMENT TABLES

## 30"x78" PRO-OPTIONS TREATMENT TABLE WITH CUSTOM CONFIGURATIONS



**PROTEAM**  
BY HAUSMANN

- Built with your training room needs in mind.
- 2" thick high-density foam upholstered tops.
- 3mm high-impact PVC edges on all doors and drawer fronts.

SKU	W	L	LAMINATE	VINYL	OPTIONS
A9075-3078-L01V19-A69A62A61	30"	78"	Natural Oak	Navy Blue	Gas Spring Backrest, Elevating Split Legs, Double Doors (6)
A9075-3078-L01V23-A62A63A61E	30"	78"	Natural Oak	Black	Gas Spring Backrest, Elevating Split Legs, Double Doors (4), Pass-Thru Modality in Middle, Emboss Logo
A9078-3078-L01V19-A69A62A61	30"	78"	Natural Oak	Navy Blue	Gas Spring Backrest, Elevating Split Legs, Double Doors (3)
A9078-3078-L01V23-A69A62A61E	30"	78"	Natural Oak	Black	Gas Spring Backrest, Elevating Split Legs, Double Doors (3), Emboss Logo

# PROTEAM™ TREATMENT TABLES

## 30"x78" H-BRACE TREATMENT TABLE WITH STORAGE SHELF



**PROTEAM**  
BY HAUSMANN

- Six (6) legs and plywood top for added strength.
- Laminate shelf with two long stretchers.
- 30" x 78" upholstered top with 2" thick high-density foam.

SKU	W	L	LAMINATE	VINYL
A9087-3078-L01V23	30"	78"	Natural Oak	Black



# PROTEAM™ TREATMENT TABLES

## 30"x78" H-BRACE SPLIT LEG TREATMENT TABLE WITH STORAGE SHELF



**PROTEAM**  
BY HAUSMANN

- Extremely versatile and functional split leg design.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Six (6) legs and full plywood subtop for added strength.

SKU	W	L	LAMINATE	VINYL
A9088-3078-L01V23	30"	78"	Natural Oak	Black

## 27"x48" H-BRACE TAPING TABLE

**PROTEAM**  
BY HAUSMANN

- 36" is convenient height for taping.
- H-brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.
- High-pressure laminate legs, apron and stretchers.



SKU	W	L	LAMINATE	VINYL
A9093-2748-L01V23	27"	48"	Natural Oak	Black

# PROTEAM™ TREATMENT TABLES

## 27"x48" H-BRACE TAPING TABLE WITH STORAGE SHELF AND OPEN CABINET STORAGE



**PROTEAM**  
BY HAUSMANN

- 36" is convenient height for taping.
- End facing open storage cabinet.
- H-brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.

SKU	W	L	LAMINATE	VINYL
A9097-2748-L01V23	27"	48"	Natural Oak	Black

## 30"x78" ELECTRIC HI-LO TREATMENT TABLE WITH SPLIT LEG AND SLIDING CABINET DOORS

**PROTEAM**  
BY HAUSMANN

- Durable motorized operation with hands-free foot control.
- Electro-hydraulic cylinders are built into four (4) black powder coated steel legs with levelers.
- Fully enclosed oak laminate cabinet with sliding laminate doors.



SKU	W	L	LAMINATE	VINYL
A9718-3078-L00V23	30"	78"	Natural Oak	Black

# WOOD TREATMENT TABLES

## H-BRACE WOOD TREATMENT TABLE



- Natural wood finish is standard.
- H-Brace tables have triple-bolted, solid hardwood legs.
- Extra-firm 2" thick high-density foam tops with rounded corners for comfort.

SKU	W	L	WOOD	VINYL	OPTIONS
ATT-2472-L00V23-A20	24"	72"	Natural Oak	Black	Upholstered Shelf
ATT-2772-L00V23-A20A69A27	24"	72"	Natural Oak	Black	Upholstered Shelf, Gas Spring Backrest, Paper Dispenser & Cutter Combo

## 30"x72" WOOD TREATMENT TABLE WITH STORAGE CABINET

- Natural wood finish is standard.
- Cabinet Tables are quality hardwood plywood and veneer construction.
- Four (4) door cabinet with center divider and one adjustable shelf on each side.



SKU	W	L	WOOD	VINYL	OPTIONS
ATT-3072-L00V20-A69	30"	72"	Natural Oak	Dove Gray	Gas Spring Backrest

# HTT TABLES

## PREMIUM H-BRACE TREATMENT TABLE - “HTT”



- 2 ½" tabletop featuring rounded corners and firm foam padding.
- High-quality Formica® laminate in oak & maple.
  - Laminate withstands heavy cleaning and disinfecting with non-degradable finish.
  - Laminate resistant to scratches and impacts under normal use conditions.
- Durable heavyweight Pro-Form vinyl upholstery – also withstands heavy cleaning and disinfecting.

SKU	W	L	LAMINATE	VINYL	OPTIONS
HTT-3072-L01V23-A09A03	30"	72"	Natural Oak	Black	U-bend Backrest, Treatment Table Shelf
HTT-3078-L01V23-A09A07A03	30"	78"	Natural Oak	Black	U-bend Backrest, Treatment Table Shelf, Nose Hole Cut-out
HTT-3078-L01V23-A09A03	30"	78"	Natural Oak	Black	U-bend Backrest, Treatment Table Shelf
HTT-3078-L01V19-A03A27	30"	78"	Natural Oak	Navy Blue	Paper Dispenser & Cutter Combo, Treatment Table Shelf



130 Union Street  
Northvale NJ 07647

1-888-HAUSMANN

[info@hausmann.com](mailto:info@hausmann.com) | [www.hausmann.com](http://www.hausmann.com)