



MAMMOTH

BY HAUSMANN

METAL TABLES
2-Section • 3-Section • Tilt



3

Three ways Mammoths are different:

WIDTH: Wider tops to accommodate a variety of patient sizes

WEIGHT: Durable design provides higher weight capacity from 440 lbs to 550 lbs*

WARRANTY: 3-years on the switches and upholstery, 5-years on the motor and 15-years on the steel frame providing long-term reliability

*Weight capacity is dependent on table model



— THE MAMMOTH DIFFERENCE

The Mammoth Collection vs The Competition

SPECIFICATIONS	COMPANY X	MAMMOTH 2 & 3	COMPANY Y	MAMMOTH T
Width	27"	30"	30"	30"
Length	76"	76"	81 1/2"	78"
Min Height	18"	18"	25 3/4"	26"
Max Height	37"	37"	35 1/8"	46"
Switch	Dual Bar	Dual Foot Controls	Hand & Foot Control	Hand Control
Weight Capacity	400 lbs	550 lbs	350 lbs	440 lbs
Foam	2"	2"	2"	2"
Motor Warranty	18 Months	5 Years	2 Years	5 Years

MAMMOTH 3



3-Section Hi-Lo Treatment Table - 7703



3-SECTION TABLE

MAMMOTH 2



2-Section Hi-Lo Examination Table - 7702



2-SECTION TABLE

MAMMOTH T

Hi-Lo Tilt Table - 7706



COLORS AND FEATURES

METAL COLOR

Cool Gray

VINYL COLORS

V19-Navy Blue

V26-Royal Blue

V20-Dove Gray

V27-Purple

V21-Forest Green

V28-Regimental Blue

V22-Deep Wine

V29-Bright Blue

V23-Black

V30-Brown

V24-American Beauty

V31-Orange

V25-Claret

V32-Pencil Gray

Key Features & Options

- Width: 30-inch wide upholstered top
- Weight: Durable design provides higher weight capacity from 440 lbs to 550 lbs*
- Warranty: 3-years on the switches and upholstery, 5-years on the motor and 15-years on the steel frame providing long-term reliability
- Tables take up less operational space with no end-to-end movement
- Heavy duty independently lockable casters
- Convenient and easy set-up
- Durable heavyweight Pro-Form vinyl upholstery - withstands heavy cleaning and disinfecting (see manual for more information)
- 2.5" tabletop featuring rounded corners and firm foam padding
- Available in both Standard and Custom SKU configurations



TILT TABLE



130 Union Street
Northvale, NJ 07647

1-888-HAUSMANN
info@hausmann.com
www.hausmann.com