

MAMMOTH 1

1-Section Hi-Lo Treatment Table - 7701



MAMMOTH
BY HAUSMANN

Introducing the new Hausmann Mammoth 1 Treatment Table. It has a durable steel frame and features high-quality Pro-Form vinyl. It is proudly assembled and upholstered in the USA. A new addition to our table portfolio with the same quality and long-lasting durability you expect from Hausmann.



Upholstered top 30" wide and 76" long

Extra width to accommodate more patients

Highly durable and long-lasting steel frame

15-year warranty on frame

No end-to-end movement

Perfect for tight spaces

550-pound weight capacity

Lifts a wider range of patients

7701* Custom SKU shown



= WIDTH & WEIGHT & WARRANTY

MAMMOTH 1

1-Section Hi-Lo Treatment Table - 7701



MAMMOTH
BY HAUSMANN

METAL COLOR

Cool Gray

VINYL COLORS

V19-Navy Blue

V20-Dove Gray

V21-Forest Green

V22-Deep Wine

V23-Black

V24-American Beauty

V25-Claret

V26-Royal Blue

V27-Purple

V28-Regimental Blue

V29-Bright Blue

V30-Brown


V31-Orange

V32-Pencil Gray

Key Features & Options

- Width: 30-inch wide upholstered top
- Weight: 550-pound weight capacity
- Warranty: 3-years on the switches and upholstery, 5-years on the motor and 15-years on the steel frame providing long-term reliability
- Two (2) foot switches for placement around the table
- Tables take up less operational space with no end-to-end movement
- Heavy-duty independently lockable casters
- 2.5" tabletop featuring rounded corners and firm foam padding
- Durable heavyweight Pro-Form vinyl upholstery - withstands heavy cleaning and disinfecting (see manual for more information)
- Convenient and easy set-up
- Optional nose cutout

Specifications



SKU	W	L	METAL	VINYL
7701*	30"	76"	Cool Gray	Custom